

# Readme: Please follow this template when writing notes

2022年8月16日 12:27

Please include:

- Date
- Name(s)/Initial(s) of the operator(s)
- Name(s) of the protein(s)/gene(s) etc. you're working on
- Name(s) of the kit(s)
  - Please specify the brand! e.g., TIANGEN Plasmid Miniprep
- Procedure (no details needed if a protocol is present)
- Any data/picture you've obtained, e.g., the picture of the gel
- Anything you want to add, e.g., 'there seems something wrong with the band, it's too dim/the size is incorrect', etc.

in your notes, so that other people can easily follow what you've done today.

Thanks!

RYY

Initials

RYY = Ru Yunyi

SC = Shen Chang

WYX/Lisa = Wu Yixuan

IJ/JJY/Irene = Jin Jingyun

# 7.28 SDS-PAGE Practice

Monday, July 25, 2022 2:59 PM

Operators: WYX/SC/JQZ/YY/RYY

Note taken by: RYY

Sample Used: Cell pellets with some His-tagged protein from last year's iGEM project, 1.5mL tubes noted 'H3.1沉淀' 1&2

Procedure:

1. Preparing Sample Lysis
  - Add 100uL 裂解缓冲液 (from GST Purification Kit) into each tube
  - Add 2uL lysosome (50mg/mL from ? kit) into each tube, mix well by pipetting and vortex
  - Let it stand on ice for 30 minutes
  - Centrifuge the samples at 15000g, 4°C
  - Aspirate the supernatant, discard the pellet
  - Add 25uL sample + 25uL loading buffer to make a total of 50uL
  - Put the samples on dry bath (90°C, 5mins)
2. SDS-PAGE Electrophoresis
  - Samples were loaded in this order: ladder ('多彩预染' \_\_\_\_), sample1, (N/A), ladder, sample2, ladder
  - Run at 120V for ~1-2 hrs
  - Stain with Coomassie Brilliant Blue (protocol by Yang)
  - Shake for 5 minutes
  - Wash with MiliQ Water for 5 minutes (change water in between)
  - Imaging using the ChemiDoc Imaging system

## 8.10 LB Plate making + bacterial culture

Tuesday, August 9, 2022 1:43 PM

LB plate making+ LB liquid medium making

LB + Kanamycin plate making (For pET30a transformation from last year to test the competence transformation efficiency)

Total media: 1L LB liquid, 200mL LB agar, 200mL LB+Amp (final conc: 20ug/mL)

Plate made by Lisa, Irene, Yunyi, Chang

## 8.15 Plasmid extraction

Monday, August 15, 2022 11:07 AM

OD600 of E. coli cultivated since August 14<sup>th</sup> 8pm

Measured at 10:15am, diluted \*10

	20IpaD	IpaD-HA	JPS-G3
LB+amp	0.147	0.149	0.145
LB	0.104	0.162	0.133

Measured at 1pm, diluted\*10

	20IpaD	IpaD-HA	JPS-G3
LB+amp	0.161	0.170	0.162
LB	0.115	0.167	0.135

OD measured by Chang

Plasmid extracted by Lisa, Chang, Yuxin (TIANGEN Plasmid Miniprep)

## 8.16 Plasmid Enzyme Digestion

Tuesday, August 16, 2022 13:44

### Double Digestion

#### [Protocol](#)

Operator: SC, RYY

Nanodrop result of extracted plasmid and the reaction system

Name	μg/μl	Applied amount in the reaction system	Enzyme(1μl)	Enzyme(1μl)	Buffer (2μl)
20IpaD_digested	0.085	11.8μl	XhoI	SacI	rcutsmart
IpaD-HA_digested	0.092	10.8μl	BamHI	XhoI	rcutsmart
JPS-G3_digested	0.119	8.4μl	XhoI	SacI	rcutsmart

Add milliQ to a reaction system of total 20 uL

\*Placed in dry bath at 37°C for 1 h for heat inactivation

Put back in -20°C for storage

## 8.22 Gel Electrophoresis & Purification

2022年8月22日 14:12

Operator: RYY


Digested plasmids 20ipaD\_digested, ipaD-HA\_digested, JPS-G3\_digested

### Gel Electrophoresis

1% agarose gel (TAE), run at 120V for 60 min

ladder	20ipaD	IpaD_HA	JPS-G3	ladder
--------	--------	---------	--------	--------

Remaining samples put back to -20°C for storage (<1uL left)



\*20ipaD missing due to running time too long  
Another double digestion session needed

### Gel Purification

#### Protocol

Beyotime kit

IpaD-HA: 124.2 mg

JPS\_G3: 113.5 mg

55°C dry bath for 5 mins

16.0xg, 1 min, 23°C

Extracted ~20uL DNA

Put in -20°C

Named as

- JPS-G3 purified 220822
- ipaD-HA purified 220822



## 8.23 Gel Electrophoresis

Tuesday, August 23, 2022 10:28

Gel electrophoresis

Operator: Chang (RYY nearby)


DNA sample: Purified IpaD-HA, Purified JPS-G3

No result

Digestion & Gel electrophoresis

Operator: Chang Shen

#	Name	μg/μl	Applied amount in the reaction system (μl)	Enzyme(1μl)	Enzyme(1μl)	Buffer (2μl)	Milli Q (μl)
1	20IpaD_digested	0.094	10.6	BamHI	SacI	rcutsmart	5.4
2	IpaD-HA digested	0.079	12.6	BamHI	XhoI	rcutsmart	3.4
3	JPS-G3_digested	0.127	7.9	XhoI	SacI	rcutsmart	10.1



Conclusion: The sample lanes were contaminated by the DNA ladder because of improper loading method, so the result was not reliable.





## 8.24 Gel Electrophoresis, Gel Purification, & Enlarge amplification of E. coli with csgA-smt gene (from XJTLU)

Wednesday, August 24, 2022 12:40

### Gel Electrophoresis

Gel made by Qizhou

Operator: Chang


Samples:

Digested Samples (9 $\mu$ l): 20IpaD, IpaD-HA, JPS-G3

Loaded Ladder (1.5  $\mu$ l \* 2 lanes): Beyotime

Loaded mixture: 10 $\mu$ l

Result:



### Gel Purification

[Protocol](#)

Beyotime kit

Operator: RYY

20IpaD: 47.7 mg


IpaD\_HA: 31.6 mg

JPG\_G3: 41.9 mg

failed, AGAIN. What's the reason????

Maybe the agarose isn't fully melted?

ladder 20IpaD IpaD-HA JPS-G3 ladder



*E. coli* (ssgA-smt, from XJTLU, No sign of bacteria on the plate after 15 h

## 8.30 Double Digestion + Gel Purification

2022年8月30日 10:29

Operator: RYY

### Double Digestion

#	Name	Concentration (µg/µl)	Applied amount in the reaction system (µl)	Enzyme-1 (1µl)	Enzyme-2 (1µl)	Buffer (5µl)	MilliQ (µl)
1	20IpaD_digested	0.139	7.2	BamHI	SacI	rcutsmart	35.8
2	IpaD-HA digested	0.113	8.8	BamHI	XhoI	rcutsmart	34.2
3	JPS-G3_digested	0.198	5.1	SacI	XhoI	rcutsmart	37.9

Incubate at 37°C for 1 hr

Incubate at 65°C for 20 mins

35uL run gel electrophoresis

### Gel Extraction

By ThermoFisher GeneJET Gel Extraction Kit

Conc:

1(20IpaD) - 0.7ng/uL;

2(IpaD\_HA) - 1.9ng/uL;

3(JPS-G3) - 1.7ng/uL

## 8.31 PCR

2022年8月31日 11:06

Operator: RYY

### PCR + Gel Electrophoresis

#### Primer pairs:

10/11 used for ipaD, 13/15 used for JPS-G3

10: CGCGGATCCATGAATATTACAACCTC (ipaD\_FW with BamHI cleavage site)

11: GGCCTCGAGTCAAGCGTAG (ipaD\_RV with XhoI cleavage site)

13: GGCCTCGAGTTATTTGTCGTCGTCATCCTTG (JPS-G3\_RV with XhoI cleavage site)

15: CGCGGATCCGGAAGTACACAAGTACAGCTAG (JPS-G3\_FW with BamHI cleavage site)

#### Protocol:


##### Reaction system:

- 0.2uL DNA template
- 1uL FP
- 1uL RP
- 12.5uL TIANGEN Taq 2xPCR Mix
- Add ddH<sub>2</sub>O to a reaction system of 25uL in total

##### Run at Bio-Rad Thermo Cycler

- 95°C for 3 min
- 98°C for 20 sec -
- 60°C for 15 sec | \*35
- 72°C for 15 sec -
- 72°C for 15 sec
- 4°C hold

##### Run at 1% agarose Gel, 69 A for ~1h



## 9.1 PCR product Purification + Double Digestion

2022年9月1日 18:15

Operator: RYY

### PCR Purification

Beyotime PCR clean-up Kit

After purification, dsDNA concentration: ipaD: 117.3ng/uL, JPS-G3: 57.7 ng/uL

### Double Digestion

#	Name	Concentration (μg/μl)	Applied amount in the reaction system (μl)	Enzyme-1 (1μl)	Enzyme-2 (1μl)	Buffer (5μl)	MilliQ (μl)
1	lpaD-HA ppcr	0.117	10	BamHI	XhoI	rcutsm art	33
2	JPS-G3 ppcr	0.057	20	BamHI	XhoI	rcutsm art	23

incubate on 37°C dry bath overnight (9.1 20:00-9.2 12:00, 16 hrs in total) + 65°C for 20 min

## 9.2 Transformation + Plate Preparation

2022年9月2日 17:18

Operator: RYY

### Transformation


pET30a from Huang Lab (KanR)

BeyoTime DH5a competent cell (stored at  $-80^{\circ}\text{C}$ )

### Protocol

To step 4, overnight culture (IB 2052)

### LB Plate Preparation



Bottle 1: 500mL MilliQ, 12.5g LB broth, 7.5 g agar

Bottle 2: 500mL MilliQ, 12.5g LB broth

Shake well

Sterilize at  $120^{\circ}\text{C}$  for 30 min

For bottle 1, add 0.4mL Sangon Kanamycin Sulfate Solution, mix well, make agar plates

Stored at  $-4^{\circ}\text{C}$  fridge

## 9.3 Strains fetched from XJTLU/Coat Plate & Overnight Culture

Saturday, September 3, 2022 19:54

Operator: RYY/JJY

### Coat Plate & Overnight Culture

Use overnight cultured pET30a to coat plates

Protocol: [Transformation](#)

1 liquid culture (5mL LB + 4uL Kan)

7 plates (1 streak plate, 1 resuspension plate, 5 serial dilution from 10X to 10<sup>-5</sup>X)

Remaining DH5a stored in -80°C in IB 2052 (labelled 'E')

---

Operator: SC

### Strains fetched from XJTLU

pGEX-4T-1-CsgA-MBP3 E. coli LB liquid culture, stored in 15ml centrifuge tube at 4 °C

pGEX-4T-1-CsgA-shMT E. coli Glycerol stock, stored at -20°C at 3080

The CsgA-shMT strain was inoculated in LB medium for recovery at 8pm.



pGEX4T-1-cs  
gA-shMT



pGEX-4T-1-C  
sgA-MBP3

About pGEX-4T-csgA-SmtA strain which was failed to culture:

The first LB+amp liquid medium inoculated with the csgA-smtA strain (for recovery) was considered to be problematic by yy as the medium didn't show that the strain was growing and the result of the plate was not normal. The medium was then stored at 4°C before it was taken out again and cultured on. I placed the medium in the incubator again on August 30.

On August 31 I checked it and found the liquid to be opaque, and placed it in 4°C.




## 9.4 Monoclonal Culture + Plasmid Extraction

2022年9月4日 22:10

Operator: RYY

### Monoclonal Culture

9:59AM single colony from overnight culture plates found, transplanted to 2\* 5mL liquid culture (LB + Kan)




### Plasmid Extraction

9:30PM start (culture time: ~12 hrs)

[Protocol: TIANGEN Plasmid MiniPrep Kit](#)

Bacterial Solution used:

- 2\* 5mL Liquid culture made today (labelled **pET30a 1 & 2**),
- polyclonal culture made yesterday (labelled **pET30a - non-MC**)



75uL/tube, stored at -20

# 9.5 Streak Plate

2022年9月5日 10:29


Operator: SC



## Plate streaking

Plate: LB+amp (made on June 2)

Strains	Source
pGEX-4T-1-Csga-MBP3	LB liquid culture from XJTLU
pGEX-4T-1-Csga-shMT	LB liquid, cultured before for 36 h
pGEX-4T-1-Csga-smtA	LB+amp liquid, from the tube considered to be problematic

Picking colony (at 12pm)





## Plate streaking

Plate: LB+amp (made by Qizhou on August 26)(Time: September 6, 1am)

Strains	Source
pGEX-4T-1-Csga-MBP3	LB liquid culture from XJTLU
pGEX-4T-1-Csga-shMT	LB liquid, cultured before for 36 h
pGEX-4T-1-Csga-smtA	LB+amp liquid, from the tube considered to be problematic

\*As suggested by yy, I used more reliable plates and streaked them again

## 9.6 Double Digestion, Gel Extraction (Failed)

2022年9月6日 22:09

Operator: RYY

### Double Digestion

#	Name	Concentration (µg/µl)	Applied amount in the reaction system (µl)	Enzyme-1 (2µl)	Enzyme-2 (2µl)	Buffer (10µl)	MilliQ (µl)
1	pET30a-1	0.07	30	BamHI	XhoI	rcutsm art	56

[protocol](#)

### Gel Extraction

Failed again..... I won't use Beyotime Kit anymore


### Double Digestion

Same as above, I will try ligation tmr

## 9.7 Gel Extraction

2022年9月7日 16:40

Operator: RYY



### Gel Extraction

pET30a: 10.8ng/uL

## 9.8 Transformation

2022年9月9日 9:01

Operator: RYY

pUC19 + JPS-G3

pUC19 + ipaD

Transformed into E. coli DH5a, overnight culture (8:40p.m.)



2022年9月9日 11:50

What we have:

1. **Target Sequences**
  - a. ipaD\_HA sequence, purified PCR product (F: BamHI/R: XhoI)
  - b. JPS-G3\_FLAG sequence, purified PCR product (F: BamHI/R: XhoI)
2. **Empty Vectors**
  - a. pET30a empty vector (KanR, LacI/Lac Operon)
  - b. pET30a empty vector digested, purified (F: BamHI/R: XhoI) (KanR, LacI/Lac Operon)
3. **Sequences in plasmids**
  - a. ipaD\_HA + pUC19 (insertion site unknown, synthesized by Sangon)
  - b. JPS-G3 + pUC19 (insertion site unknown, synthesized by Sangon)
  - c. pGEX-4T-1-CsgA-MBP3 (from XJTLU)
  - d. pGEX-4T-1-CsgA-shMT (from XJTLU)
4. **E.coli strains**
  - a. DH5a competent cell (T7 polymerase-)
  - b. BL21(DE3) competent cell (T7 polymerase+)
5. **Primers**

## Goal 1: Express & Purify proteins ipaD\_HA & JPS-G3\_FLAG, Test Affinity

- ☒ **Ligate**<sub>(1)</sub> (1.a + 2.b) (1.b + 2.b) to make pET30a\_ipaD\_HA and pET30a\_JPS\_G3  
☒ Use **colony PCR**<sub>(2)</sub> (two plasmids transformed into DH5a competent cells) to test whether the plasmid is intact, then  
☐ **Transform**<sub>(3)</sub> useful plasmids to BL21(DE3) competent cells, and  
☐ **Induce**<sub>(4)</sub> the expression of target proteins via IPTG, then  
☐ **Purify**<sub>(5)</sub> proteins, and  
☐ protein **affinity analysis via ELISA**<sub>(6)</sub>

## Goal 2: Bind CsgA to ipaD\_HA & JPS-G3 (i.e., Surface Display)

- ☐ Overlap extension PCR (or called gene SOE):


CsgA+Linker with ipaD\_HA = *Bam*HI\_CsgA+Linker+ipaD\_HA\_ *Xho*I


- 1st round PCR:
  - Tube 1: pGEX-4T-1-CsgA-shMT + (OvF\_csgA-shMT + OvR1\_csgA-shMT), Tm = 63
  - Tube 2: pET30a\_ipaD\_HA + (OvF\_ipaD + 一开始订的Primer11), Tm = 51
- 2nd round PCR:
  - Tube 1 product + Tube 2 product + (OvF\_csgA-shMT + OvR\_ipaD), Tm = 61

Expected Final Product:



csgA-linker-i  
paD-3xHA





分区 Group E 的第 26 页

## 9.9 Ligation + Transformation

2022年9月9日 12:51

Operator: RYY

### Ligation

<b>pET30a</b>	<b>6uL (~50ng, 0.014pmol)</b>
<b>JPS-G3</b>	<b>0.3uL (~14ng, 0.042pmol)</b>
10X Ligation Buffer	<b>2μl</b>
双蒸水或Milli-Q水	<b>11.2uL</b>
T4 DNA ligase	<b>0.5uL</b>
总体积	20μl

<b>pET30a</b>	<b>6uL (~50ng, 0.014pmol)</b>
<b>ipaD</b>	<b>0.3uL (~30ng, 0.042pmol)</b>
10X Ligation Buffer	<b>2μl</b>
双蒸水或Milli-Q水	<b>11.2uL</b>
T4 DNA ligase	<b>0.5uL</b>
总体积	20μl

JPS-G3: 57ng/uL, ipaD: 117ng/uL

Incubate at 25 from 15:30 to 17:00 (1.5 hr), then stored at -20

Biomass calculator: <https://www.promega.com.cn/resources/tools/biomath/>

### Transformation

Transformed into DH5a competent cells, overnight culture

No results

## 9.12 Purification of Plasmids and Ligation

Monday, September 12, 2022 21:56

Operator: RYY

### Purification of Plasmids and Insertion Sequences

Double digestion of plasmid again, purified by both GeneJET gel extraction kit (pET30a-2) and Beyotime DNA purification kit (pET30a-1; proved to be better);  
used Beyotime DNA purification kit to purify insertion sequences again

#### dsDNA Concentration

Conc.	ng/uL
ipaD (~1000bp)	41.5
JPS-G3 (~600bp)	35.4
pET30a-1	7.9
pET30a-2	16.9

### Overnight Ligation

	Empty vector	Insert sequence	T4	Ligation Buffer	DNAse-free water	Result
i-1	1: 6.5uL	ipaD: 1uL	0.5uL	2uL	10uL	++
J-1	1: 6.5uL	JPS-G3: 1uL	0.5uL	2uL	10uL	++
i-2	2: 3uL	ipaD: 1uL	0.5uL	2uL	13.5uL	-
J-2	2: 3uL	JPS-G3: 1uL	0.5uL	2uL	13.5uL	+

Plasmid: 0.015 pmol

Insert sequence: 0.05 pmol

Overnight culture at 16°C

## 9.13 Transformation


Tuesday, September 13, 2022 20:42

Operator: RYY

J-1, J-2, I-1, I-2 transformed into DH5a, overnight culture

Result:







## 9.14 Colony PCR + Overnight Culture; Picking colonies + overnight culture for strains from XJTLU


2022年9月14日 16:22

Operator: RYY

### Colony PCR

4 colonies were picked for each gene, run PCR + Gel electrophoresis

All positive




Operator: Chang


Two colonies were picked for strain *csgA-smtA*

Four colonies were picked for strain *csgA-smht* (one by RYY)

Plate:

*csgA-smtA*





csgA-smht





## 9.15 Plasmid extraction, enzyme digestion and Gel electrophoresis for csgA; ipaD & JPS-G3 recombinant plasmids sent for sequencing

Thursday, September 15, 2022 13:51

Operator: Chang Shen

### Plasmid Extraction

Source	Plasmid concentration (ng/ $\mu$ l)
smht-1	66.7
smht-2	95.4
smht-3	98.2
smht-4(RYY)	69.1
smtA-1	99.0
smtA-2	41.8

### Enzyme digestion

Applied Enzyme: BamHI & XhoI

Predicted length:

csgA-smtA	656bp
csgA-smht	656bp


Reaction system:


Plasmid	Plasmid volume ( $\mu$ l)	BamHI ( $\mu$ l)	XhoI( $\mu$ l)	rCutSmart	MilliQ ( $\mu$ l)	Total ( $\mu$ l)
smht-1	15.00	1	1	5	8	30
Smht-2	10.48	1	1	5	12.52	30
Smht-3	10.18	1	1	5	12.82	30
Smht-4(RYY)	14.47	1	1	5	8.53	30
smtA-1	10.10	1	1	5	12.9	30
smtA-2	23.93	1	1	5	4.07	30

Gel electrophoresis result

Supposed length

csgA-shMT	656bp
csgA-smtA	656bp













Operator: RYY

## Plasmid Extraction & Sequencing

ipaD-1 (D1) & ipaD-3 (D2) & JPS-G3-5 (J1) & JPS-G3-6 (J2) Plasmids extracted, sent for sequencing

Name	Conc (ng/uL)	Sequencing result	.ab file (T7)	.ab file (T7t)
<b>D1</b>	60.2	Point Mutation A>G	 G220901353 _D1_T7_E03	 D1
<b>D2</b>	60.4	No problem	 G220901354 _D2_T7_F03	 D2
<b>J1</b>	58.3	No problem	 G220901355 _J1_T7_A07	 G220901355 _J1_T7t_B07
<b>J2</b>	57.7	Point mutation A>G	 G220901356 _J2_T7_C07	 G220901356 _J2_T7t_D07

# 9.19 Overlap Extension PCR

2022年9月19日 15:55

Operator: RYY

## PCR 1

Goal: Amplify csgA for

1. F/R1: ipaD (expected length: 497 bp)
2. F/R2: JPS-G3 (expected length: 498 bp)

Templates used: shmt-2/shmt-3 plasmids

### Reaction system:

- 0.5uL DNA template
- 1uL FP
- 1uL RP
- 12.5uL TIANGEN Taq 2xPCR Mix
- Add ddH2O (10uL) to a reaction system of 25uL in total

### Protocol:

- 95°C for 3 min
- 98°C for 20 sec -
- 58°C for 15 sec | \*35 <- 58 = 63-5
- 72°C for 15 sec -
- 72°C for 15 sec
- 4°C hold

Purified by Beyotime PCR product purification kit

## PCR 2

Goal: Amplify JPS-G3 (expected length: 531 bp)

Templates used: pET30a\_JPS-G3

### Reaction system:


- 0.5uL DNA template
- 1uL FP
- 1uL RP
- 12.5uL TIANGEN Taq 2xPCR Mix
- Add ddH2O (10uL) to a reaction system of 25uL in total

### Protocol:

- 95°C for 3 min
- 98°C for 20 sec -
- 54°C for 15 sec | \*35 <- 54 = 59-5
- 72°C for 15 sec -
- 72°C for 15 sec
- 4°C hold

Purified by Beyotime PCR product purification kit

csgA (expected: 497bp)/JPS-G3 (expected: 531bp), length error



## 9.20 Overlap Extension PCR

2022年9月22日 15:23

Operator: RYY

### PCR 1

Goal: Amplify ipaD (expected length: 1109bp)

Templates used: pET30a-ipaD


#### Reaction system:

- 0.4uL DNA template
- 1uL FP
- 1uL RP
- 12.5uL TIANGEN Taq 2xPCR Mix
- Add ddH<sub>2</sub>O (10uL) to a reaction system of 25uL in total

#### Protocol:

- 95°C for 3 min
- 98°C for 20 sec -
- 50°C for 15 sec | \*35
- 72°C for 15 sec -
- 72°C for 15 sec
- 4°C hold

Purified by Beyotime PCR product purification kit



### PCR 2

Goal: Overlap Extension

Templates used: product from 19-PCR1, 19-PCR2, 20-PCR1

#### Reaction system:

- 0.4uL + 0.4uL DNA templates
- 1uL FP
- 1uL RP
- 12.5uL TIANGEN Taq 2xPCR Mix
- Add ddH<sub>2</sub>O (10uL) to a reaction system of 25uL in total

**Protocol:**

- 95°C for 3 min
- 98°C for 20 sec -
- 55°C for 15 sec | \*35
- 72°C for 15 sec -
- 72°C for 15 sec
- 4°C hold